



Sutter Medizintechnik GmbH • Tullastr. 87 • 79108 Freiburg/Germany

TO WHOM IT MAY CONCERN

22. April 2016 /ps

Confirmation: Minimum and maximum time for our validated sterilization method

To whom it may concern:

We hereby confirm that the exposure time of 3 minutes mentioned in our IFU (DOC_REF 89 90 00) for the sterilization method (Fractionated vacuum) reflects the validated minimum time for a successful sterilization.

Please note: The maximum exposure time for this sterilization method (Fractionated vacuum) is 5 minutes and the maximum exposure temperature is 138°C

With best regards,

Sutter Medizintechnik

i.A. Patric Schlegel

